





Dyna Filters Pvt. Ltd.

Sterilizing Tunnel



DYNA MAKE TUNNEL



Model	Overall dimensions (mm)		
	Length	width	Height
Steri Dyn 450	3735	1015	2340
Steri Dyn 600	3960	1220	2340
Steri Dyn 900	3960	1525	2340

Tech-brief	
Governing standards	ISO 14644 (1998)
AIR velocity	90 fpm (0.45mps) \pm 20% in all zones except sterilization, where it is maintained in the range 130-160 fpm (0.65 - 0.80 mps)
Cleanliness	Fewer than 100 particles of size 0.5 or larger per cu. Ft of air sampled (or < 3,520/cuM)