

STATIC PASS BOX

TECHNICAL SPECIFICATIONS: -

MATERIAL OF CONSTRUCTION: - S.S.304/316 PVC coated matt finish OR CRCA powder coated with RAL9002 colour shade.

ACCESSORIES: -

Hour meter connected to UV light.
Door handles.

ELECTRICALS: -

8/20w, Fluorescent light with milky diffuser.
8/15w, UV light along with Hour meter.
12Volt, Electromagnet for door interlocking.
12 Volt, Buzzer.
12 Volt Indicator lamps.
16 Amps. Three-pin top along with connection cable.

SELECTION CHART: -

MODEL NO	OVERALL DIMENSIONS (MM)			WORKSPACE DIMENSIONS (MM)			ELECTRICALSUPPLY (PHASE)
	W	D	H	W	D	H	SINGLE
DFPL/SPB/2-2-2	730	690	760	610	610	610	SINGLE
DFPL/SPB/3-3-3	1035	995	1055	915	915	915	SINGLE

NOTE – ALL OTHER DIMENSIONS ARE AVAILABLE AS PER REQUIREMENT.